


Differential Gas Diffusion Porometer

Applications

Pore structure characterization of membranes containing very small pores can be performed. Delicate membranes sensitive to pressure can also be tested. Pore diameter as well as flow rates are measurable.


Principle

The flow rate of a gas is measured as a function of applied differential pressure across the sample. The gas flows through all the pores larger than the pore diameter of the gas molecule. The test is repeated using gases having different molecular diameters. The flow rates are usually very small. The sample is loaded, the system is evacuated, and the inlet chamber is maintained at a constant pressure. Time rate of pressure increase in the outlet chamber is measured. The flow rate, F , through the sample in volume at standard temperature and pressure is computed.




$$F = [V_o T_s / p_s T] (dp/dt)$$

where p is gas pressure on the outlet chamber, t is time, V_o is volume of the outlet chamber, T_s is standard temperature, T is test temperature and p_s is standard pressure. Flow rate as low as $10^{-4} \text{ cm}^3/\text{s}$ is measurable.



Change of outlet gas pressure with time for a sample



Measurable Characteristics

- Pore throat diameter
- Bubble point & Mean flow pore diameter
- Flow Distribution


For samples having very small pores, mean free path of gas is greater than pore diameter and gas flow is primarily by diffusion through pores larger than the molecular diameter of gas.

$$[F(M)^{1/2} / \Delta p] = (1/6p_s l)(2\pi RT)^{1/2} \sum D^3$$

F = flow rate at p_s and T, M = molecular weight of gas, D = pore diameter, Δp = pressure drop, l = thickness, R = gas constant
Distribution function; f

$$f = -dF'/dD$$

where F' is cumulative flow function, $[F(M)^{1/2} / \Delta p]$



Features

- Throat diameters of through pores measurable
- Very small pore diameters measurable
- Flow rates through the small pores measurable
- Flow distribution measurable
- Usable over wide pressure and temperature intervals employing a wide variety of gases
- No damage to product
- Fully automated, very little operator involvement and robust

Other Products

Average Fiber Diameter Analyzer

Bubble Point Tester

Capillary Flow Porometer

Capillary Condensation Flow Porometer

Complete Filter Cartridge Analyzer

Clamp-On Porometer

Compression Porometer

Custom Porometer

Cyclic Compression Porometer

Envelope Surface Area Analyzer

Filtration Media Analyzer

High Flow Porometer

Integrity Analyzer

In-Plane Porometer

Microflow Porometer

Nanopore Flow Porometer

QC Porometer

Diffusion Permeameter

Gas Permeameter

Liquid Permeameter

Vapor Permeameter

Water Vapor Transmission Analyzer

Liquid Extrusion Porosimeter

Mercury/Nonmercury Intrusion Porosimeter

Vacuapore

Water Intrusion Porosimeter (Aquapore)

BET Liquisorb

BET Sorptometer

Gas Pycnometer

Mercury Pycnometer

Also Available:

Testing Services

Consulting Services

Short Courses

Buy Rent Lease

Porous Materials, Inc.

20 Dutch Mill Rd, Ithaca, NY 14850 USA

Tel: (607)-257-5544 Toll Free in USA & Canada: 1-800-TALK-PMI

Fax: (607) 257-5639 Email: info@pmiapp.com WWW.PMIAPP.COM

