

The Complete Filter Cartridge Analyzer

Applications

Fully assembled filter cartridges are widely used in many industries including biotechnology, pharmaceutical, chemical, beverage, and food. Filtration efficiency in all these applications is determined by the pore structure of the complete filter cartridge. The Filter Cartridge Analyzer measures the bubble point, the mean flow pore diameter, and the pore distribution of the complete cartridge rather than a small sample of the filter media. The tester also measures the gas permeability.


Instrument

The unique instrument design easily deals with very high volume of gas flowing through large cartridges and reduces pressure drop by eliminating narrow ducts, bends and constrictions. A tank for storage of gas under pressure sufficient for the test is supplied as part of the test equipment so that standard laboratory air supply is adequate for test execution. The sample chamber holds the cartridge between a fixed head and an adjustable head. By adjusting the position of the adjustable head cartridge of any length can be accommodated. A pneumatically operated piston applies sufficient pressure to seal the edges of the cartridge.



Principle

A wetting liquid is allowed to spontaneously fill the pores of the cartridge and air pressure is increased to empty pores and permit gas flow. Measured differential gas pressure and flow rates through a cartridge in wet and dry conditions yield various pore structure characteristics.



Capability

Pore diameter: The measured differential gas pressure yields the through pore throat diameter.

$$D = 4 \gamma \cos \theta / p$$

D = pore diameter

γ = surface tension of wetting liquid

θ = contact angle of the liquid

p = differential gas pressure.



Bubble point: Computed from pressure for initiation of flow through wet sample.

Mean flow pore diameter: Computed from mean flow pressure at which wet curve and half-dry curves meet.

Pore distribution: Given in terms of distribution function, $f \cdot f = -d[(F_w/F_d) \times 100]/dD$

where F_w and F_d are wet and dry flow respectively.

Gas permeability: Computed from the gas flow rates using Darcy's law.



Features

- Completely automated
- Windows based software, simple operation, and minimal operator involvement
- Only a few minutes for test execution
- Sintered metal, woven metal, polymeric, and ceramic cartridges can be tested
- Adequate safety precautions

Other Products

Average Fiber Diameter Analyzer
 Bubble Point Tester
 Capillary Flow Porometer
 Capillary Condensation Flow Porometer
 Complete Filter Cartridge Analyzer
 Clamp-On Porometer
 Compression Porometer
 Custom Porometer
 Cyclic Compression Porometer
 Envelope Surface Area Analyzer
 Filtration Media Analyzer
 High Flow Porometer
 Integrity Analyzer

In-Plane Porometer
 Microflow Porometer
 Nanopore Flow Porometer
 QC Porometer
 Diffusion Permeameter
 Gas Permeameter
 Liquid Permeameter
 Vapor Permeameter
 Water Vapor Transmission Analyzer
 Liquid Extrusion Porosimeter
 Mercury/Nonmercury Intrusion Porosimeter
 Vacuapore
 Water Intrusion Porosimeter (Aquapore)

BET Liquisorb
 BET Sorptometer
 Gas Pycnometer
 Mercury Pycnometer

Also Available:
 Testing Services
 Consulting Services
 Short Courses

Buy Rent Lease

Porous Materials, Inc.
 20 Dutch Mill Rd, Ithaca, NY 14850 USA
 Tel: (607)-257-5544 Toll Free in USA & Canada: 1-800-TALK-PMI
 Fax: (607) 257-5639 Email: info@pmiapp.com

WWW.PMIAPP.COM

