

Membrane Manufacturing Apparatus

Some of the hollow fiber membranes made in the PMI membrane manufacturing machines were tested by Liquid-Liquid Porometry. The presence of small pores in the membrane are demonstrated by the plot of flow rate with pressure (Figure 6). The mean flow pore diameter was 0.015 microns. Most of the pores are present over a narrow size range as demonstrated by the pore distribution in Figure 7.

Figure 6.

Flow rates through dry and wet hollow fiber samples PVDF Polymer





Figure 7.

Pore Distribution of hollow fiber membranes PVDF Polymer



Specifications

Typical Properties for Hollow Fiber based on Polymer Type:

Mean pore size: 0.01-50 microns
ID, OD and Wall Thickness: Dependent on Spinnerette polymer before OD, Needle OD, Needle ID
Length: 500 ft plus

Typical Properties for Cast membrane based on Polymer Type:

Mean pore size: 0.01-100 microns
Thickness: 0.001" - 0.20"
Length: 3-4 ft
Width: 6" max

Hollow Fiber Manufacturing Machine



Flat Sheet Membrane Manufacturing Machine

Application

PMI membrane manufacturing machines can be used to create Cast flat sheet membranes and hollow fiber membranes. The machines permit adjustment of fabrication parameters so that membranes with different characteristics could be made for development, research, and many wide varieties of applications and filtration.



July 2011

20 Dutch Mill Rd, Ithaca, NY 14850, USA
Toll Free (US & Canada): 1-800-TALK-PMI (1-800-825-5764)
Phone: 607-257-5544 Fax: 607-257-5639

Email: info@pmiapp.com www.pmiapp.com

Principle

Solutions for making membranes are thoroughly mixed at desired temperature and pressure. The chemicals are shaped in to flat sheets or hollow fibers and allowed to coagulate at the desired temperature.

The Machines

Layout of the hollow fiber making machine is shown in Figure 1. One hundred PSI pressure tanks hold chemicals at the desired temperature and pressure. The chemicals are constantly stirred and pumped to the spinnerette. Flow regulators control the flow rates of the chemicals. The hollow fibers coming out of the spinnerette pass through a temperature controller coagulation tank and are washed. Figure 2 shows the layout of the flat sheet membrane making machine.

Figure 1.
Hollow Fiber Making Machine





Figure 2.
Flat sheet membrane making machine.



Typical Products

A variety of hollow fiber membranes and flat sheet membranes have been manufactured. Figure 3 shows typical hollow fiber membranes.




Figure 3.
Hollow Fiber Membranes

Pore Structure of Membranes

The pore structures of the membranes made in the PMI membrane making machines were determined by various pore structure characterization techniques. The pore structures of hydrophobic flat sheet membranes were measured by water intrusion porosimetry. Figure 4 shows variation of hydrophobic pore volume with pore diameter. Pores have diameter from about, 50 to 0.01 microns. The pore distribution in Figure 5 shows that pores making maximum contribution to the pore volume are about - 0.04 microns in diameter. The pore volume at 20 micron is probably due to pore mouths.

Figure 4.
Pore Volume and Diameter of Flat Sheet Membrane PVDF Polymer




Figure 5.
Pore Distribution of Flat Sheet Membrane PVDF Polymer

